



ADULTER-SKREEN controls are manufactured to qualitatively monitor the performance of laboratory urine adulterant testing procedures. *Please read the entire package insert before using the controls. Each Kitbox contains one 20ml bottle of levels A & B.*

**INTENDED USE**

The ADULTER-SKREEN controls are intended to validate the performance of laboratory urinary adulterant test methods. They should be treated as any "unknown" specimen while following the specific protocol of the assay being used.

*This product is intended to be used by health care professionals as an integral part of good laboratory practices.*

**SUMMARY AND EXPLANATION**

The ADULTER-SKREEN product line of controls is a synthetic matrix that have been stabilized to ensure that the product will be viable until the date of expiration. Subsequent testing of the control material is performed qualitatively to ascertain that abnormal controls test abnormal and normal controls test normal for the claimed life of the product.

**PROCEDURE**

1. Allow controls to come to room temperature followed by **gentle swirling or inversion** before use.
2. Transfer an appropriate aliquot of control urine as required by the test device or screening method.

**DESCRIPTION**

ADULTER-Screen controls are manufactured in accordance with the following specifications:

Constituent	CONTROL A	CONTROL B
Creatinine	Negative	>90 mg/dL (High Normal)
Nitrite	Negative	Positive
Glutaraldehyde	Positive	Negative
pH	≤4.0 (Low Abnormal)	>9.0 (High Abnormal)
Specific Gravity	≤1.005 (Low Abnormal)	≥1.025 (High Abnormal)
Oxidant	Negative	Positive (PCC) Positive (H <sub>2</sub> O <sub>2</sub> )

**EXPECTED RESULTS**

The abnormal ADULTER-SKREEN constituents must test abnormal and the normal ADULTER-SKREEN constituents must test normal.

**STORAGE & STABILITY**

**Unopened:**

The controls are stable until the expiration date when stored refrigerated at (2-8°C) and protected from light.

**After Opening:**

The controls are stable for 60 days or until the expiration date, whichever comes first, when stored refrigerated (2-8°C).

**PRECAUTIONS**

**For in vitro diagnostic use only**

Please read the entire package insert before using the ADULTER-SKREEN control urines. Please use the same safety precautions you would use for processing any "unknown" sample. Protect product from exposure to direct sunlight. *Do not use beyond the expiration date.*

**LIMITATIONS OF PROCEDURE**

This control is meant to be used to validate the performance of qualitative urine adulterant screening methods.

Consult test manufacturer's instructions when using this product; changes in reagents, sample requirement, or methodology may affect test results. Due to the variable nature of colorimetric test pads on adulterant devices and test strips by various manufacturers, certain constituents color may not match exactly to color intensity indicated on chart.

*This product is not meant to be used as a standard or calibrator.*



**Biochemical Diagnostics, Inc.**  
**180 Heartland Blvd.**  
**Edgewood, NY 11717 USA**

Phone: (631) 595-9200  
Email: support@biochemicaldiagnostics.com  
Website: www.kovaintl.com

For additional information on our other products please contact us or refer to our website.

SYMBOL LEGEND	
	Consult Instructions for Use
	Temperature Limits
	In Vitro Diagnostic Medical Device
	Batch Code
	Product Catalog Number
	Unique Device Identifier
	Manufacturers Identification
	Country of Manufacture
	Use by Date
	Caution, Consult Accompanying Documents
	For Prescription Use Only