

DETECTABUSE® STAT-SKREEN 23 AND 23H LIQUID URINE CONTROL					
ABBREVIATIONS	CONSTITUENTS	NEGATIVE	CUTOFF	STAT-SKREEN 23 POSITIVE	STAT-SKREEN 23H POSITIVE
THC	Delta-9-THC-COOH	None Detected	50	100	100
COC	Benzoyllecgonine	None Detected	300	600	600
PCP	PCP (Phencyclidine)	None Detected	25	50	50
MOR/OPI	Morphine	None Detected	300 / 2000	600	4000
AMP	d-Amphetamine	None Detected	1000	2000	2000
MET	d-Methamphetamine	None Detected	1000	2000	2000
BAR	Secobarbital	None Detected	300	600	600
BUT/BAR	Butalbital	None Detected	300	600	600
BZO	Oxazepam	None Detected	300	600	600
METH	Methadone	None Detected	300	600	600
PPX	Propoxyphene	None Detected	300	600	600
TCA	Nortriptyline	None Detected	1000	2000	2000
MDMA	MDMA*	None Detected	500	1000	1000
OXY	Oxycodone	None Detected	100	200	200
BUP	Buprenorphine	None Detected	10	20	20
BUPG	Buprenorphine Glucuronide	None Detected	10	20	20
FEN	Norfentanyl	None Detected	Varies	50	50
TRAM	Tramadol	None Detected	Varies	400	400
6MAM	6-MAM	None Detected	10	20	20
ETG	ETG	None Detected	Varies	1000	1000
EDDP	EDDP	None Detected	300	600	600
K2	JWH-018 5- Pentanoic acid	None Detected	50	100	100

*3,4 Methylenedioxy methyl-amphetamine